

Anti-Human/Mouse/Rat DAP Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09774

3th Edition

Description: This gene encodes a basic, proline-rich, 15-kD protein. The protein acts as a positive mediator of programmed cell death that is induced by interferon-gamma.

Antigen: Synthetic peptide of human DAP

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: