

Anti-Human/Mouse/Rat DDAH2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09779

3th Edition

Description: This gene belongs to the dimethylarginine dimethylaminohydrolase (DDAH) gene family. The encoded enzyme plays a role in nitric oxide generation by regulating cellular concentrations of methylarginines, which in turn inhibit nitric oxide synthase activity.

Antigen: Synthetic peptide of human DDAH2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Detected in heart, placenta, lung, liver, skeletal muscle, kidney and pancreas, and at very low levels in brain.