

Anti-Human/Mouse/Rat KDM8 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09823

3th Edition

Description: This gene likely encodes a histone lysine demethylase. Studies of a similar protein in mouse indicate a potential role for this protein as a tumor suppressor. Alternatively spliced transcript variants have been described.

Antigen: Recombinant protein of human KDM8

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Weakly expressed in most cells. Highly expressed in breast cancer cells.