

Anti-Human/Mouse JMJD7 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09827

3th Edition

Description: This gene encodes a highly conserved protein with a JmjC domain, which are part of the cupin metalloenzyme superfamily. JmjC proteins may function as 2-oxoglutarate-Fe(II)-dependent dioxygenases. Most tissues also express read-through transcripts from this gene into the downstream phospholipase A2, group IVB (cytosolic) gene, some of which Human, Mouse encode fusion proteins combining the N-terminus of this protein with the phospholipase A2, group IVB protein.

Antigen: Recombinant protein of human JMJD7

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: