

Anti-Human STK17A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09835

3th Edition

Description: This gene is a member of the DAP kinase-related apoptosis-inducing protein kinase family and encodes an autophosphorylated nuclear protein with a protein kinase domain. The protein has apoptosis-inducing activity.

Antigen: Synthetic peptide of human STK17A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed in placenta. Lower levels in heart, lung, skeletal muscle, kidney and pancreas.