

Anti-Human/Mouse/Rat KCNK1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09839

3th Edition

Description: This gene encodes one of the members of the superfamily of potassium channel proteins containing two pore-forming P domains. The product of this gene has not been shown to be a functional channel, however, it may require other non-pore-forming proteins for activity.

Antigen: Recombinant protein of human KCNK1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed with high levels in heart and brain and lower levels in placenta, lung, liver and kidney.