

**Anti-Human/Mouse/Rat KLHL2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09856**

---

3th Edition

**Description:**KLHL2 (kelch-like protein 2), also known as actin-binding protein Mayven, is a 593 amino acid actin-binding protein. Highly expressed in brain, KLHL2 is thought to play a role in the organization of the actin cytoskeleton in brain cells.

**Antigen:**Recombinant protein of human KLHL2

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Predominantly expressed in brain.