

Anti-Human/Mouse/Rat FEM1B Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09919

3th Edition

Description: This gene encodes an ankyrin repeat protein that belongs to the death receptor-associated family of proteins and plays a role in mediating apoptosis. The encoded protein is also thought to function in the replication stress-induced checkpoint signaling pathway via interaction with checkpoint kinase 1

Antigen: Synthetic peptide of human FEM1B

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed. Highly expressed in testis. Weakly expressed in other tissues.