

**Anti-Human/Mouse/Rat LZTFL1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09924**

---

3th Edition

**Description:** This gene encodes a ubiquitously expressed protein that localizes to the cytoplasm. This protein interacts with Bardet-Biedl Syndrome (BBS) proteins and, through its interaction with BBS protein complexes, regulates protein trafficking to the ciliary membrane. Nonsense mutations in this gene cause a form of Bardet-Biedl Syndrome; a ciliopathy characterized in part by polydactyly, obesity, cognitive impairment, hypogonadism, and kidney failure. This gene may also function as a tumor suppressor; possibly by interacting with E-cadherin and the actin cytoskeleton and thereby regulating the transition of epithelial cells to mesenchymal cells. Alternative splicing of this gene results in multiple transcript variants.

**Antigen:** Recombinant protein of human LZTFL1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in prostate, ovary, stomach, pancreas, esophagus, breast, liver, bladder, kidney, thyroid, colon and lung (at protein level). Down-regulated in multiple primary tumors (at protein level). Detected in testis, heart, skeletal muscle, thymus, spleen, small intestine, and peripheral blood leukocytes.