

Anti-Human MAGEB10 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09932

3th Edition

Description: This gene encodes a member of the B subfamily of the melanoma associated antigen protein family. The encoded protein is specifically expressed in testis and tumor cells.

Antigen: Recombinant protein of human MAGEB10

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: