

Anti-Human/Mouse TUBGCP2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA09977

3th Edition

Description:Gamma-tubulin complex component 2 is a protein that in humans is encoded by theTUBGCP2 gene. It is part of the gamma tubulin complex, which required formicrotubule nucleation at the centrosome. Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome. Gamma-tubulin complex is composed of gamma-tubulin, TUBGCP2, TUBGCP3, TUBGCP4, TUBGCP5 and TUBGCP6.

Antigen:Synthetic peptide of human TUBGCP2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitously expressed.