

**Anti-Human/Mouse NDUFAF2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09995**

---

3th Edition

**Description:** NADH:ubiquinone oxidoreductase (complex I) catalyzes the transfer of electrons from NADH to ubiquinone (coenzyme Q) in the first step of the mitochondrial respiratory chain, resulting in the translocation of protons across the inner mitochondrial membrane. This gene encodes a complex I assembly factor. Mutations in this gene cause progressive encephalopathy resulting from mitochondrial complex I deficiency.

**Antigen:** Recombinant protein of human NDUFAF2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Highly expressed in ESCC cells. Also expressed in heart, skeletal muscle, liver, and in fibroblasts.