

**Anti-Human A1BG Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10047**

---

3th Edition

**Description:** Alpha-1-B glycoprotein is a protein in humans that is encoded by the A1BG gene. The protein encoded by this gene is a plasma glycoprotein of unknown function. The protein shows sequence similarity to the variable regions of some immunoglobulin supergene family member proteins.

**Antigen:** Synthetic peptide of human A1BG

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Plasma.