

Anti-Human TRPA1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10093

3th Edition

Description:The structure of the protein encoded by this gene is highly related to both the protein ankyrin and transmembrane proteins. The specific function of this protein has not yet been determined; however, studies indicate the function may involve a role in signal transduction and growth control.

Antigen:Synthetic peptide of human TRPA1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed at very low level. Expressed at very low level in human fibroblasts and at a moderate level in liposarcoma cells.