

**Anti-Human/Mouse/Rat PITPNM1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10117**

---

3th Edition

**Description:**PITPNM1 belongs to a family of membrane-associated phosphatidylinositol transfer domain-containing proteins that share homology with the Drosophila retinal degeneration B (rdgB) protein.

**Antigen:**Recombinant protein of human PITPNM1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Ubiquitous.