

**Anti-Human/Mouse/Rat NIT1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10119**

---

3th Edition

**Description:** This gene encodes a member of the nitrilase protein family with homology to bacterial and plant nitrilases, enzymes that cleave nitriles and organic amides to the corresponding carboxylic acids plus ammonia. Multiple transcript variants encoding different isoforms have been found for this gene.

**Antigen:** Recombinant protein of human NIT1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Detected in heart, brain, placenta, liver, skeletal muscle, kidney and pancreas.