

Anti-Human/Mouse/Rat ID3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10138

3th Edition

Description:The protein encoded by this gene is a helix-loop-helix (HLH) protein that can form heterodimers with other HLH proteins. However, the encoded protein lacks a basic DNA-binding domain and therefore inhibits the DNA binding of any HLH protein with which it interacts.

Antigen:Synthetic peptide of human ID3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed abundantly in lung, kidney and adrenal gland, but not in adult brain.