

Anti-Human OCIAD2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10163

3th Edition

Description:OCIAD2 (OCIA domain containing 2) is a 154 amino acid protein that localizes to endosomes and contains one OCIA domain. OCIAD2 exists as multiple alternative spliced isoforms and is encoded by a gene which localizes to human chromosome 4. Representing approximately 6% of the human genome, chromosome 4 contains nearly 900 genes, one of which is the Huntingtin gene, which is found to encode an expanded glutamine tract in cases of Huntington's disease.

Antigen:Recombinant protein of human OCIAD2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: