

Anti-Human/Mouse/Rat SLC32A1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10193

3th Edition

Description: Vesicular inhibitory amino acid transporter is a protein that in humans is encoded by the SLC32A1 gene. The protein encoded by this gene is an integral membrane protein involved in gamma-aminobutyric acid (GABA) and glycine uptake into synaptic vesicles. The encoded protein is a member of amino acid/polyamine transporter family II.

Antigen: Synthetic peptide of human SLC32A1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Retina. Expressed throughout the horizontal cells or more specifically at the terminals.