

**Anti-Human PDIA2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10277**

---

3th Edition

**Description:**Protein disulfide isomerases (EC 5.3.4.1), such as PDIP, are endoplasmic reticulum (ER) resident proteins that catalyze protein folding and thiol-disulfide interchange reactions

**Antigen:**Synthetic peptide of human PDIA2

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Highly expressed in pancreas (at protein level).