

**Anti-Human/Mouse/Rat RMDN3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10298**

---

3th Edition

**Description:**Involved in cellular calcium homeostasis regulation. May participate in differentiation and apoptosis of keratinocytes. Overexpression induces apoptosis. Present at high level in epidermis and seminiferous epithelium: while basal cells in the epidermis and spermatogonia show no perceptible amount, keratinocytes of suprabasal layers and differentiating first-order spermatocytes up to spermatids exhibit high expression.

**Antigen:**Recombinant protein of human RMDN3

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Present at high level in epidermis and seminiferous epithelium: while basal cells in the epidermis and spermatogonia show no perceptible amount, keratinocytes of suprabasal layers and differentiating first-order spermatocytes up to spermatids exhibit high expression. In skeletal muscle, its presence is restricted to fibers of the fast twitch type. In surface epithelia containing ciliated cells, it is associated with the microtubular structures responsible for ciliary movement. Also present in specific structures of the central nervous system such as neurons of the hippocampal region, ganglion cells of the autonomic nervous system, and axons of the peripheral nervous system (at protein level). Widely expressed.