

Anti-Human/Mouse/Rat PTPN5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10311

3th Edition

Description:Protein tyrosine phosphatase non-receptor type 5 is an enzyme that in humans is encoded by the gene. May regulate the activity of several effector molecules involved in synaptic plasticity and neuronal cell survival, including MAPKs, Src family kinases and NMDA receptors.

Antigen:Recombinant protein of human PTPN5

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: