

Anti-Human/Mouse/Rat RAB9A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10332

3th Edition

Description: Ras-related protein Rab-9A is a protein that in humans is encoded by the RAB9A gene. Rab 9A is a 201 amino acid protein that localizes to the cytoplasmic side of the cell membrane, as well as to the membrane of the Golgi apparatus and the ER, and is involved in the transport of proteins between endosomes and the trans Golgi network.

Antigen: Recombinant protein of human RAB9A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: