



**Anti-Human/Mouse/Rat RBBP9 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10365**

---

3th Edition

**Description:** The protein encoded by this gene is a retinoblastoma binding protein that may play a role in the regulation of cell proliferation and differentiation. Two alternatively spliced transcript variants of this gene with identical predicted protein products have been reported, one of which is a nonsense-mediated decay candidate.

**Antigen:** Recombinant protein of human RBBP9

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Widely expressed. Expressed at higher levels in tumor tissues than in normal tissues.