

**Anti-Human/Mouse/Rat NDRG3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10397**

---

3th Edition

**Description:**The N-myc downstream regulated gene (NDRG) family is comprised of four members, NDRG1 (also designated Drg1, RTP, rit42, Cap43 and Ndr1), NDRG2, NDRG3 and NDRG4, which share 57-65% homology. The gene encoding NDRG3 maps to human chromosome 20q11 and is predominantly expressed in testis, prostate and ovary, which suggests it may play a role in spermatogenesis.

**Antigen:**Synthetic peptide of human NDRG3

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Ubiquitous. Highly expressed in brain.