

Anti-Human RNASET2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10400

3th Edition

Description: This ribonuclease gene is a novel member of the Rh/T2/S-glycoprotein class of extracellular ribonucleases. It is a single copy gene that maps to 6q27, a region associated with human malignancies and chromosomal rearrangement.

Antigen: Recombinant protein of human RNASET2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitous. Higher expression levels observed in the temporal lobe and fetal brain.