

Anti-Human/Mouse MVD Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10465

3th Edition

Description:The enzyme mevalonate pyrophosphate decarboxylase catalyzes the conversion of mevalonate pyrophosphate into isopentenyl pyrophosphate in one of the early steps in cholesterol biosynthesis. It decarboxylates and dehydrates its substrate while hydrolyzing ATP.

Antigen:Synthetic peptide of human MVD

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in heart, skeletal muscle, lung, liver, brain, pancreas, kidney and placenta.