

**Anti-Human/Mouse ADRB3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10530**

---

3th Edition

**Description:**The protein encoded by this gene belongs to the family of beta adrenergic receptors, which Human, Mouse mediate catecholamine-induced activation of adenylate cyclase through the action of G proteins. This receptor is located mainly in the adipose tissue and is involved in the regulation of lipolysis and thermogenesis.

**Antigen:**Synthetic peptide of human ADRB3

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed mainly in adipose tissues.