

**Anti-Human/Mouse SLC25A26 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10534**

---

3th Edition

**Description:** Mitochondrial carriers, including SLC25A26, are a family of transport proteins found mostly in the inner membranes of mitochondria. They shuttle metabolites and cofactors through the mitochondrial membrane

**Antigen:** Recombinant protein of human SLC25A26

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Widely expressed. Highly expressed in testis, with moderate expression in brain, heart, kidney, lung, skeletal muscle, pancreas, small intestine and liver, and low expression in spleen.