

**Anti-Human SMYD4 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10555**

---

3th Edition

**Description:**The SMYD4 antibody from Proteintech is a rabbit polyclonal antibody to a fusion protein of human SMYD4. This antibody recognizes human, mouse, rat antigen. The SMYD4 antibody has been validated for the following applications: ELISA, WB, IHC analysis.

**Antigen:**Recombinant protein of human SMYD4

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**