

Anti-Human/Mouse/Rat MYCN Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10559

3th Edition

Description: This gene is a member of the MYC family and encodes a protein with a basic helix-loop-helix (bHLH) domain. This protein is located in the nucleus and must dimerize with another bHLH protein in order to bind DNA. Amplification of this gene is associated with a variety of tumors, most notably neuroblastomas.

Antigen: Synthetic peptide of human MYCN

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in the neuronal cells of the cerebrum, neuroblastomas and thyroid tumors (at protein level).1
Publication