

**Anti-Human/Mouse/Rat SNX15 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10569**

---

3th Edition

**Description:** This gene encodes a member of the sorting nexin family. Members of this family contain a phosphoinositide binding domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. Overexpression of this gene results in a decrease in the processing of insulin and hepatocyte growth factor receptors to their mature subunits. This decrease is caused by the mislocalization of furin, the endoprotease responsible for cleavage of insulin and hepatocyte growth factor receptors.

**Antigen:** Recombinant protein of human SNX15

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Widely expressed.