

**Anti-Human/Mouse/Rat SDCBP Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10598**

---

3th Edition

**Description:**The protein encoded by this gene was initially identified as a molecule linking syndecan-mediated signaling to the cytoskeleton. The syntenin protein contains tandemly repeated PDZ domains that bind the cytoplasmic, C-terminal domains of a variety of transmembrane proteins. This protein may also affect cytoskeletal-membrane organization, cell adhesion, protein trafficking, and the activation of transcription factors.

**Antigen:**Recombinant protein of human SDCBP

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Widely expressed. Expressed in fetal kidney, liver, lung and brain. In adult highest expression in heart and placenta.