

Anti-Human/Mouse/Rat STX7 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10649

3th Edition

Description:The Syntaxin 7 antibody from Proteintech is a rabbit polyclonal antibody to a recombinant protein of human Syntaxin 7. This antibody recognizes human,mouse,rat antigen. The Syntaxin 7 antibody has been validated for the following applications: ELISA, WB, IHC, IP analysis.

Antigen:Recombinant protein of human STX7

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Highest expression is found in placenta followed by heart, skeletal muscle, kidney and brain. Low expression is found in pancreas, lung and liver.