

**Anti-Human/Mouse SOCS7 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10749**

---

3th Edition

**Description:**The protein as regulator signaling cascades probably through protein ubiquitination and/or sequestration. It plays an important role in insulin signaling and glucose homeostasis through IRS1 ubiquitination and subsequent proteasomal degradation. Inhibits also prolactin, growth hormone and leptin signaling by preventing STAT3 and STAT5 activation, sequestering them in the cytoplasm and reducing their binding to DNA.

**Antigen:**Synthetic peptide of human SOCS7

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in brain and leukocytes. Also in fetal lung fibroblasts and fetal brain.