

Anti-Human/Mouse/Rat RAB5A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10761

3th Edition

Description: Ras-related protein Rab-5A is a protein that in humans is encoded by the RAB5A gene. It is rate-limiting for homotypic early endosome fusion. It required for the fusion of plasma membranes and early endosomes, as well as contributes to the regulation of filopodia extension.

Antigen: Synthetic peptide of human RAB5A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: