

Anti-Human/Mouse RAB5C Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10764

3th Edition

Description:The protein is a member of the Rab protein family.Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment .The funtion is protein transport,Probably involved in vesicular traffic.

Antigen:Synthetic peptide of human RAB5C

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: