



**Anti-Human BMP3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10803**

---

3th Edition

**Description:**BMP3 belongs to the transforming growth factor-beta (TGFB) superfamily. Bone morphogenic protein, also known as osteogenin, induces bone formation.

**Antigen:**Synthetic peptide of human BMP3

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in adult and fetal cartilage.