

Anti-Human SAMD4A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10899

3th Edition

Description: Sterile alpha motifs (SAMs) in proteins such as SAMD4A are part of an RNA-binding domain that functions as a posttranscriptional regulator by binding to an RNA sequence motif known as the Smaug recognition element, which was named after the *Drosophila* Smaug protein.

Antigen: Synthetic peptide of human SAMD4A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: