



**Anti-Human/Mouse SAMD4B Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10900**

---

3th Edition

**Description:** SMGB; hSmaug2

**Antigen:** Synthetic peptide of human SAMD4B

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Widely expressed in embryonic and adult tissues.1 Publication