

**Anti-Human/Mouse/Rat CWC27 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA10927**

---

3th Edition

**Description:**CWC27 belongs to the cyclophilin-type PPlase family and contains 1 PPlase cyclophilin-type domain. There are two isoforms. PPlases accelerate the folding of proteins.

**Antigen:**Synthetic peptide of human CWC27

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**