

Anti-Human/Mouse/Rat CTGF Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA10996

3th Edition

Description: The protein encoded by this gene is a mitogen that is secreted by vascular endothelial cells. The encoded protein plays a role in chondrocyte proliferation and differentiation, cell adhesion in many cell types, and is related to platelet-derived growth factor. Certain polymorphisms in this gene have been linked with a higher incidence of systemic sclerosis.

Antigen: Synthetic peptide of human CTGF

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in bone marrow and thymic cells. Also expressed one of two Wilms tumors tested.