

Anti-Human SLC25A38 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11028

3th Edition

Description: This gene is a member of the mitochondrial carrier family. The encoded protein is required during erythropoiesis and is important for the biosynthesis of heme. Mutations in this gene are the cause of autosomal congenital sideroblastic anemia.

Antigen: Synthetic peptide of human SLC25A38

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Preferentially expressed in erythroid cells.