

**Anti-Human SLC40A1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11049**

---

3th Edition

**Description:**The protein encoded by this gene is a cell membrane protein that may be involved in iron export from duodenal epithelial cells. Defects in this gene are a cause of hemochromatosis type 4 (HFE4).

**Antigen:**Synthetic peptide of human SLC40A1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in placenta, intestine, muscle and spleen.