

Anti-Human/Mouse/Rat SLIT1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11087

3th Edition

Description: Slit homolog 1 protein is a protein that in humans is encoded by the SLIT1 gene. Thought to act as molecular guidance cue in cellular migration, and function appears to be mediated by interaction with roundabout homolog receptors. During neural development involved in axonal navigation at the ventral midline of the neural tube and projection of axons to different regions. By similarity, SLIT1 and SLIT2 together seem to be essential for midline guidance in the forebrain by acting as repulsive signal preventing inappropriate midline crossing by axons projecting from the olfactory bulb.

Antigen: Synthetic peptide of human SLIT1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Predominantly expressed in adult forebrain. Expressed in fetal brain, lung and kidney.