

Anti-Human/Mouse SMC4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11113

3th Edition

Description:Members of the structural maintenance of chromosomes, or SMC, family (e.g., SMC1A; MIM 300040) are critical for mitotic chromosome condensation in frogs and for DNA repair in mammals.

Antigen:Synthetic peptide of human SMC4

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Widely expressed. Higher expression in testis, colon, thymus.