

**Anti-Human HINT2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11146**

---

3th Edition

**Description:** Histidine triad proteins, such as HINT2, are nucleotide hydrolases and transferases that act on the alpha-phosphate of ribonucleotides. Hydrolase probably involved in steroid biosynthesis. May play a role in apoptosis. Has adenosine phosphoramidase activity. High expression in liver and pancreas. Expression is significantly down-regulated in hepatocellular carcinoma (HCC) patients.

**Antigen:** Synthetic peptide of human HINT2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** High expression in liver and pancreas. Expression is significantly down-regulated in hepatocellular carcinoma (HCC) patients.