

**Anti-Human/Mouse/Rat TMSB4X Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11224**

---

3th Edition

**Description:** This gene encodes an actin sequestering protein which plays a role in regulation of actin polymerization. The protein is also involved in cell proliferation, migration, and differentiation. This gene escapes X inactivation and has a homolog on chromosome Y.

**Antigen:** Synthetic peptide of human TMSB4X

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in several hemopoietic cell lines and lymphoid malignant cells. Decreased levels in myeloma cells.