

Anti-Human/Mouse TRIM45 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11262

3th Edition

Description: This gene encodes a member of the tripartite motif family. The encoded protein may function as a transcriptional repressor of the mitogen-activated protein kinase pathway. Alternatively spliced transcript variants have been described.

Antigen: Synthetic peptide of human TRIM45

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in skeletal muscle, brain, heart and pancreas.