

Anti-Human/Mouse/Rat FOXJ3 Polyclonal Antibody
Polyclonal Antibody

Cat.NO.: PA11268

3th Edition

Description: forkhead/winged helix transcription factor, Foxj3, which was expressed in embryonic and adult skeletal muscle. The 100-kb, 13-exon mouse Foxj3 gene on chromosome 4 encodes a 623 amino acid (aa) protein from an mRNA of at least 4.8 kb (Human FOXJ3: Chr 1, 627 aa, 5.3-kb mRNA). During the stages of mouse development investigated Foxj3 is expressed in neuroectoderm, in neural crest, and in many structures derived from neural crest cells, such as facioacoustic, trigeminal, and dorsal root ganglia. Stripes of expression appear at E10.5 in the location of myotomes and expand ventrally in a pattern similar to the developing body wall musculature. Developing limbs have a complex pattern of Foxj3 expression that at E12.5 colocalizes with the condensed mesenchyme of the skeletal primordia.

Antigen: Synthetic peptide of human FOXJ3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: